



Axona

Product Catalogue

(Rev. 1.1)
15/10/2015

*Axona Ltd
Unit 4U St. Albans Enterprise Centre
Long Spring, Porters Wood
St. Albans, Herts, AL3 6EN
United Kingdom*

<http://www.axona.com>
support@axona.com

Table of Contents

Digital Data Acquisition Systems	3
Additional Components	4
Software	5
Tethers	6
Headstages	8
Adaptors	10
Microdrives	11
Re-usable Standard Drives	11
Disposable Custom Drives	12
Parts & Wiring Components	13

*****PLEASE NOTE THAT ALL IMAGES CONTAINED IN THIS CATALOGUE ARE REPRESENTATIONS ONLY. CONTACT AXONA FOR MORE INFORMATION ABOUT YOUR SPECIFIC ORDER*****





DIGITAL DATA ACQUISITION SYSTEMS






Axona's core product is its multichannel data acquisition system (dacq), designed for recording single neurons from multiple indwelling electrodes in awake, behaving animals. It has an inbuilt digital oscilloscope which reduces the size and cost of the system.




Dacq systems interface with a number of inbuilt behavioural data collection systems:

- Position tracker allows tracking of LEDs on the animals, in two or three dimensions.
- Digital I/O system allows event-logging and apparatus control (eg., food pellet delivery, or monitoring of lever presses)
- Stimulator allows control of a stimulus isolator (ie., for LTP or similar kinds of experiments)

Our data acquisition systems come complete with a versatile control program which combines features for pre-recording screening of brain signals (including a ten-channel digital oscilloscope display). This role of this software is primarily to configure the recording and apparatus control systems, collect electrophysiological and/or behavioural data and store it on disk for off-line analysis.





Product Code	Description	Notes
<p><i>DacqUSB/16</i></p> 	16-channel recording system (includes hardware and software)	Can later be upgraded to more channels
<p><i>DacqUSB/32</i></p> 	32-channel recording system (includes hardware and software)	Can later be upgraded to more channels
<p><i>DacqUSB/64</i></p> 	64-channel recording system (includes hardware and software)	Can later be upgraded to more channels
<p><i>DacqUSB/128</i></p> 	128-channel recording system (includes hardware and software)	-
<i>DacqUSB/UP</i>	Recording system upgrade	-
<i>SYS/ACC-1</i>	Digital accelerometer preamp module	-
<i>SYS/ACC-2</i>	digital accelerometer sensor board and cable	-



<p>VID/TRACK</p> 	<p>10.4" Tracker monitor</p>	<p>-</p>
<p>VID/CAM</p> 	<p>Camera</p>	<p>-</p>
<p>VID/WEBCAM</p>	<p>Webcam</p>	<p>For use with USB</p>
<p>VID/LENS</p> 	<p>CS-mount varifocal lens</p>	<p>-</p>
<p>VID/MULTI</p>	<p>Two-camera multiplexor</p>	<p>-</p>
<p>VID/DIGI</p>	<p>USB video digitiser for DacqTrack</p>	<p>-</p>
<p>SYS/TEST</p> 	<p>Test signal generator</p>	<p>For testing headstages, cables, system settings, etc</p>
<p>ADAPT/TEST</p>	<p>Mill-Max test signal generator adaptor</p>	<p>For use with SYS/TEST</p>
<p>SYS/STIM</p> 	<p>Stimulus isolator</p>	<p>Required to deliver external stimulus to subject via implant.</p>
<p>DacqUSB/Isol1</p> 	<p>DacqUSB isolator</p>	<p>Optional intermediary used to isolate the preamp where ground loops are an issue (e.g. foot shocks, etc).</p>
<p>SYS/REP</p> 	<p>Input Replicator</p>	<p>"Copies" 16 channels to two inputs. This lets you record units (e.g. on ch1-16) and EEG (on ch17-32) from the same electrodes.</p>

<p>COMP/TINT01</p> 	<p>TINT cluster-cutting analysis software license</p>	<p>Per license fee.</p>
<p>COMP/TINT02</p> 	<p>Supplemental TINT software license</p>	<p>Per license fee. Can only be purchased after an initial COMP/TINT01 purchase</p>
<p>COMP/TRACK</p> 	<p>DacqTrack software license</p>	<p>Per license fee.</p>

TETHERS

Fine Wire tethers are used to connect a headstage to an Axona recording system. Tethers of 8 to 32 channels are available. They can be purchased in a standard 1 or 3 meter length to minimise noise-pickup ahead of the preamplifier, but custom lengths and extension tethers can be supplied.

Product Code	Description	Input/Output
<p><i>CAB/8-3M</i></p> 	8-channel headstage cable with connectors (3m)	Input: 44-socket high density D connector (female) Output: 12-pin MillMax connector (male)
<p><i>CAB/8-1M</i></p> 	8-channel extension cable with connectors (1m)	Input: 44-socket high density D connector (female) Output: 12-pin MillMax connector (male)
<p><i>CAB/16-3M</i></p> 	16-channel headstage cable with connectors (3m)	Input: 44-socket high density D connector (female) Output: 20-pin MillMax connector (male)
<p><i>CAB/16-1M</i></p> 	16-channel extension cable with connectors (1m)	Input: 44-socket high density D connector (female) Output: 20-pin MillMax connector (male)

<p><i>CAB/32-3M</i></p> 	<p>32-channel headstage cable with connectors (3m)</p>	<p>Input: 44-socket high density D connector (female) Output: Two 20-pin MillMax connectors (male)</p>
<p><i>CAB/32-1M</i></p> 	<p>32-channel extension cable with connectors (1m)</p>	<p>Input: 44-socket high density D connector (female) Output: Two 20-pin MillMax connectors (male)</p>
<p><i>CAB/CUST</i></p>	<p>Custom specification cable</p>	<p>-</p>

HEADSTAGES

Axona headstage amplifiers are miniature, unity gain buffer amplifier circuit with AC-coupled inputs which are intended to be located directly at the recording site, and isolate the preparation from the cable connecting it to the recording system.

Although there are two sizes of boom available on directional headstages, the smaller version is highly recommended as it appears to work just as well in normal ceiling height environments, and is less inclined to get caught on things and damaged. For similar reasons it is also recommended that the boom be oriented parallel to the connector (two attachment points for the boom with a bit of spread between them is stronger).




Specifications:

Unity gain: 1.00000

Power for op-amps: +/- 5V

Product Code	Description	Dimensions (mm) L x W x H	Weight (g)
HS/A8-DIR50 HS/A8-DIR82 	8-channel directional Axona headstage	23 x 11 x 16 (small boom – 50, large boom – 82)	3.46
HS/MM8-DIR50 HS/MM8-DIR82 	8-channel directional MillMax headstage	27 x 17 x 3 (small boom – 50, large boom – 82)	3.04
HS/A16-DIR50 HS/A16-DIR82 	16-channel directional Axona headstage	20 x 20 x 15 (small boom – 50, large boom – 82)	5.12
HS/MM16-DIR50 HS/MM16-DIR82 	16 (to 32)-channel directional MillMax headstage (LEDs may be oriented along either axis)	24 x 24 x 3 (to 6) (small boom – 50, large boom – 82)	5.94

<p><i>HS/A8</i></p> 	<p>8-channel Axona connector headstage (no LED)</p>	<p>23 x 11 x 16</p>	<p>1.87</p>
<p><i>HS/A8-LED</i></p> 	<p>8-channel Axona connector headstage (1 LED)</p>	<p>23 x 11 x 16</p>	<p>2.1</p>
<p><i>HS/A16</i></p> 	<p>16-channel Axona connector headstage (no LED)</p>	<p>20 x 20 x 15</p>	<p>2.55</p>
<p><i>HS/A16-LED</i></p> 	<p>16-channel Axona connector headstage (1 LED)</p>	<p>20 x 20 x 15 (flexible LED – 10)</p>	<p>2.80</p>
<p><i>HS/MM16</i></p> 	<p>16-channel MillMax connector headstage (no LED)</p>	<p>20 x 24 x 3</p>	<p>2.29</p>
<p><i>HS/MM16-LED</i></p> 	<p>16-channel MillMax connector headstage (1 LED)</p>	<p>20 x 24 x 3 (flexible LED – 10)</p>	<p>2.29</p>

<p>HS/OM16</p> 	<p>16-channel Omnetics headstage with MillMax connector</p>	<p>17 x 4 x 28</p>	<p>1.98</p>
<p>ADPT/OM-MM-16</p>  <p>(dual 16ch version shown)</p>	<p>16-channel Omnetics-MillMax adapter (male or female)</p>	<p>27 x 17 x 11 (with pins)</p>	<p>2.67</p>
<p>ADPT/OM-MM-32</p> 	<p>32-channel Omnetics-MillMax adapter (male or female)</p>	<p>27 x 17 x 11 (with pins)</p>	<p>2.74</p>
<p>ADPT/OM44-MM-32</p>	<p>32ch Omnetics(44pin)-Mill-Max adapter (male or female)</p>	<p>-</p>	<p>-</p>
<p>ADPT/NANOZ</p>	<p>32-channel Omnetics-MillMax adaptor (Poor Lady to nanoZ)</p>	<p>-</p>	<p>-</p>

MICRODRIVES




We offer two main categories of chronic implants: re-usable standard and disposable custom drives. Re-usable drives come with 8 and 16 channels with an easy-fit Axona connector and 32 channels with a MillMax connector. These are attached to a sturdy steel frame and include a fixed ground wire. The tetrodes are all driven simultaneously with one screw, and the travel distance for the is 6mm.

Custom drives all have MillMax connectors and come with 8 to 64 channels. Each tetrode can be driven individually, although if you would rather drive them in groups, we can modify the design for you. Travel distance for custom drives is 4-5mm; for a fee (~£10), the drive bodies can be milled to accommodate larger screws and up to 10mm driving distance.


With the custom drives you can have single- or multi-site implants within a single drive, optogenetic fiber ports, etc. The exit points of the tetrodes can be situated in an number of places to best accommodate your implant site position.

We do not currently offer Omnetics connector drives; please see the FAQ on our website for more information about your Omnetics options with us.

****Below are a few examples of the kinds of drives we offer, but please contact us first so we can create the drive that best suits your needs!****

Product Code	Description	Dimensions (mm) L x W x H	Weight (g)
<i>DR/08</i> 	8-channel re-usable microdrive (Axona connector)	23 x 6 x 21	1.797 (wired)
<i>DR/16</i> 	16-channel re-usable microdrive (Axona connector)	23 x 9 x 25	2.34 (wired)
<i>DR/32</i> 	32-channel re-usable microdrive (MillMax connector)	22 x 14 x 22	4.07 (wired)












<p><i>DR/01</i></p> 	<p>Bare re-usable frame</p>	<p>9 x 6 x 21</p>	<p>0.71</p>
<p><i>Contact us for options</i></p>	<p>8-channel disposable custom drive (MillMax connector)</p>	<p>9 x 8 x 13.5</p>	<p>0.29 (unwired)</p>
<p><i>Contact us for options</i></p>	<p>16-channel disposable custom drive (MillMax connector)</p>	<p>12 x 9 x 17.5</p>	<p>1.385 (wired)</p>
<p><i>Contact us for options</i></p>	<p>32-channel disposable custom drive (MillMax connector)</p>	<p>12 x 14 x 17</p>	<p>2.34 (wired)</p>
<p><i>Contact us for options</i></p>	<p>64-channel disposable custom drive (MillMax connector)</p>	<p>32 x 23 x 24</p>	<p>6.46 (unwired)</p>

PARTS & WIRING COMPONENTS

When ordering drives, headstages, adaptors, etc. you will be advised as to which additional parts you need to complete your drive or equipment. If you need to order additional parts, you can use the list below to determine which ones

Product Code	Description	Notes
CON/17F	17-pin headstage connector (female)	-
CON/17M	17-pin drive connector (male)	Pre-wired
CON/9F	9-pin headstage connector (female)	-
CON/9M	9-pin drive connector (male)	Pre-wired
CON/MMM	Mill-Max male connector strip (50x2)	-
CON/MMF	Mill-Max female connector strip (50x2)	-
CON/25D	25-pin D-connector	For cable ports on preamp
CON/44D	44-pin D-connector	For I/O port on system unit
CON/D-CASE	Plastic shell case for D-connector	-
CON/OM-16M	16-channel male Omnetics connector	-
CON/OM-16F	16-channel female Omnetics connector	-
CON/OM-32M	32-channel male Omnetics connector	NPD-36-VV-GS
CON/OM-32F	32-channel female Omnetics connector	NSD-36-DD-GS

CANN/11	Drive cannulae 11mm	For use with custom drives
CANN/15	Drive cannulae 15mm	For use with custom drives
CANN/21G-12 	Shielding cannula, pack of 10	21G x 12mm For use with re-usable drives
CANN/17G-7 	Guide cannula, pack of 10	17G x 7mm For use with re-usable drives
PART/PIN01	17-pin (16ch) pin-protector	-
PART/PIN02	9-pin (8ch) pin-protector	-
PART/PIN03 	Drive gold pins	For use with custom drives
PART/PIN04 	Insect pins (pack of 100)	For use with custom drives
PART/PIN05	Drive metal dowels (pack of 100)	For use with 64ch custom drives
PART/A-PIN	Axona microdrive pins	-
PART/A-SOCK	Axona microdrive sockets	-
PART/M1.0X6 	Drive screws, M1.0 x 6	For use with custom drives (shuttles)
PART/M1.2X12 	Drive screws, M1.2 x 12	For use with custom drives (binding)
WIRE/12T 	Tungsten wire	12.5um (0.0005"), quote per 100'

<p>WIRE/17PI</p> 	<p>Platinum-Iridium wire</p>	<p>17um (0.00067"), quote per 100'</p>
<p>WIRE/25PI</p> 	<p>Platinum-Iridium wire</p>	<p>25um (0.001"), quote per 100'</p>